

The key feature of the Eye/ Face Wash is the ability of the cushion-flow atomizer heads to deliver a soft, feather-like flow effect to the injured ocular area.

RECEPTOR BOWL

High visibility yellow extra-large Atomizer Heads Cicolac Receptor Bowl resistant to damage from alkalis, salt solutions, oils, most acids.

ATOMISER HEADS

Twin cushion-flow Cicolac Atomizer Heads in high visibility Yellow/Black with flip-flop dust covers that protect heads that automatically release with water pressure

ACTIVATION

Yellow Cicolac Push Plate/ Foot Pedal operated Stay-open Brass Ball Valve with Stainless Steel Ball and Spindle.

STREAM CONTROL

Steady water flow under varying pressures assured with automatic pressure compensators for 30- 90 PSI.

SUPPLY

1/2" BSP

Recommended pressure: 30 PSI/ 2 kg/cm2

DRAIN

1 1/4" BSP

NET SHIPPING WEIGHT

LSE 1 : 8 kg

LSE 1A : 2 kg

LSE 1B : 3 kg

LSE 2 : 9 kg



**LSE 1 GI: EYE/ FACE WASH
- PUSH LEVER OPERATED**

**LSE 2 GI: EYE/ FACE WASH
- FOOT PEDAL OPERATED**



**LSE 1A GI : EYE/ FACE WASH AND BODY SPRAY
WITH 6' HOSE, SQUEEZE VALVE - WALL BRACKETED**



**LSE 1B GI: EYE/ FACE WASH PUSH LEVER
OPERATED - WALL MOUNTED**